



1. Tiles axes:

- DAPI images from your scan corresponding to the pattern and position chosen in 2. and 3.

2. Scan pattern selection boxes:

- Select the pattern corresponding to how your scan images were acquired
- See diagram below for examples

3. Change tiles in view:

- Select “show new positions” to show randomly selected tiles
- Use text box to manually specify a position to view

4. Select controls points:

- Launches MATLAB cpselect tool.
- Select matching pixels between rows and columns to guide stitching
- Exit out of cpselect window to save control points
- You can continue adding control points for multiple tiles in scan.
- Final values are averaged when closing the GUI.

5. Preview stitch:

- Use the stored control points to stitch the tiles in view.